



PelySenSe

Compact and low-noise quartz-enhanced photoacoustic sensor for sub-ppm ethylene detection in atmosphere

P. Patimisco^{a,b}

A. Sampaolo^a, M. Giglio^a, F. Sgobba^a, H. Rossmadl^c, V. Mackowiak^c, B. Gross^c, A. Cable^d, F.K. Tittel^b, and V. Spagnolo^{a,b}

^a PolySense Lab - Dipartimento Interateneo di Fisica, Politecnico and University of Bari, Italy; ^bRice University, Department of Electrical and Computer Engineering,, USA;

^cThorlabs GmbH, Dachau, Germany ^dThorlabs, Inc., Newton, USA






Outline

- Introduction to QEPAS
- Basics on QTFs operating at the 1° overtone mode
- QEPAS sensor for ethylene detection
- Results and performance
- Conclusions

Quartz-Enhanced PAS Physics and Basic Operation





Q-factor

P. Patimisco et al., "Recent advances in quartz enhanced photoacoustic sensing", Applied Physics Review, in press, 2018

P. Patimisco et al., "Quartz-enhanced photoacoustic spectrophones exploiting custom tuning forks: a review", Advances in Physics X 2, 169-187, 2016

Quartz tuning fork



Resonance frequency

Euler-Bernoulli beam theory


$$EI_y \frac{d^4X}{dz^4} + \rho S \frac{d^2X}{dt^2} = 0$$


Resonance Frequency, in-plane flexural modes

$$f_n = \frac{\pi \mathbf{T}}{8\sqrt{12}\mathbf{L^2}} \sqrt{\frac{E}{\rho}} n^2$$



 $f_{1st \ overtone} \sim 6.2 \cdot f_{fundamental}$





P. Patimisco, et al., "A quartz enhanced photo-acoustic gas sensor based on a custom tuning fork and a THz QCL," Analyst 139, 2079, 2014.

A. Sampaolo, et al., "Quartz-enhanced photoacoustic spectroscopy exploiting tuning fork overtone modes," Appl. Phys. Lett., 107, 231102, 2015

P. Patimisco et al., "Analysis of the electro-elastic properties of custom QTFs for optoacoustic gas sensing" Sens. Actuators B Chem. 227, 539, 2016

Quartz tuning fork Why overtone mode?



Euler-Bernoulli beam theory

$$X_n(z) = A((\cosh(k_n z) - \cos(k_n z)) - (\frac{\cosh(k_n L) + \cos(k_n L)}{\sinh(k_n L) + \sin(k_n L)})(\sinh(k_n z) - \sin(k_n z)))$$

Fundamental mode Overtone mode 1.0 1.0 0.5 0.5 QEPAS Signal (mV) 0.0 QEPAS Signal (mV) 0.3 0.1 12 10 12 14 Distance from the QTF top (mm) Distance from the QTF top (mm)


Quartz tuning fork Quality factor



$$\frac{1}{Q(P)} = \frac{1}{Q_{air}(P)} + \frac{1}{Q_{sup}}$$

Air Damping losses

$$Q_{air} = \frac{4\pi\rho T w^2 f_n}{3\pi\mu w + \frac{3}{4}\pi w^2 \sqrt{4\pi\rho_{air}\mu f_n}}$$



Support losses

$$Q_{sup} = A_n \frac{L^3}{T^3} \qquad \frac{A_{fund}}{A_{overt}} \sim 12$$

$$\frac{A_{fund}}{A_{overt}} \sim 12$$


Straightforward approach for reducing support losses is designing QTFs having prongs with large length-to-thickness aspect ratio.

Quartz- tuning fork Design and test



Optical test


H₂O line @1296.49 cm⁻¹




Ethylene detection Motivations

PolySenSe

- Colorless, odorless and flammable hydrocarbon.
- Produced by the **petrochemical industry** as a result of steam cracking
- In the **food industry**, measures quantity and the ripeness state of the fruit
- Human breath analysis (severity of oxidative stress and metabolic disturbances)




Where I can detect ethylene?



PolySenSe

Experimental Setup




Preliminary investigation




Pressure: 120 Torr

Flow: 20.8 sccm

Modulation:



Preliminary investigation





Pressure: 120 Torr Flow: 20.8 sccm

Selected absorption line

Wavenumber: 966.38 cm⁻¹ Linestrength: 2.2 10⁻²⁰ cm/mol

Laser setpoints:

Temperature: 15°C Current: 442.1 mA

Modulation:

Amplitude: 23 mV Frequency: 10751.59 Hz

Ramp:


Amplitude: 100 mV Frequency: 5 mHz

Integration time: 100 ms

PolySenSe

Micro-resonator tubes






| QTF frequency | prong spacing (ps) (mm) | ID (mm) | $\frac{ID}{ps}$ |
|------------------|----------------------------|------------|-----------------|
| 32.7 kHz | 0.3 | 0.6 | 2 |
| 7.5 kHz | 0.8 | 1.3 | 1.6 |
| 4.2 kHz | 1.0 | 1.5 | 1.5 |



 $ID = 0.85 \, mm$

Tube length




$$\frac{\lambda}{4} \le L \le \frac{\lambda}{2}$$



$$\frac{\lambda}{4} \le L \le \frac{\lambda}{2} \qquad \Longrightarrow \quad L = \frac{\lambda}{4} = 4 \text{mm}$$

Results





Pressure: 120 Torr

Flow: 20.8 sccm

Selected absorption line

Wavenumber: 966.38 cm⁻¹ Linestrength: 2.2 10⁻²⁰ cm/mol

Laser setpoints:

Temperature: 15°C Current: 442.1 mA

Modulation:

Amplitude: 23 mV Frequency: 10751.59 Hz


Ramp:


Amplitude: 100 mV Frequency: 5 mHz

Integration time: 100 ms

QEPAS Sensor Calibration







Linear response
Slope: 0.85 mV/ppm



Signal enhancement X 27

PolySenSe

Allan Deviation Analysis


@100 ms Integration time:

1σ-Noise: 0.086 mV Peak signal: 85.84 mV



Signal to Noise Ratio (SNR): 1000

Minimum detection limit (MDL): 100 ppb



@10 s Integration time:

1σ-Noise: 0.026 mV



MDL: 30 ppb

QEPAS Sensor Let's compact it!





Conclusions



- High-Q QTF operating at overtone mode
- QEPAS Sensor for ethylene detection
- Detection limit of 30 ppb at 10 s integration time

Future Perspectives

- New generation of QTF
- Single-tube micro-resonator system